

Ordered by: _____ Phone # (_____) _____

BILLING: P.O. Number _____ Townsend Account # _____

Bill To: _____

Ship To: _____

Address: _____

Address: _____

City: _____

City: _____

State: _____ Zip Code: _____ Country: _____

State: _____ Zip Code: _____ Country: _____

Phone: (_____) _____ Fax: (_____) _____

Phone: (_____) _____ Fax: (_____) _____

Shipping Preference: Ground 2-Day P.M. 2-Day A.M. Next Day P.M. Next Day A.M.

(If no preference is indicated, this order will be shipped 2 Day P.M.) Note: We do not ship products directly to patients.

Patient's Last Name: _____

Patient's First Name: _____

Compartment

Medial Lateral

If you select lateral OA, please read these instructions and select the lateral OA hinge position.

Optimum unloading for Lateral OA patients is achieved when the hinge is placed on the medial side of the knee. However, if the patient is extremely knock-knee and a medial hinge is problematic, request the Lateral Hinge position.

Medial Hinge Position Lateral Hinge Position

Leg: Left Right

Size

X-Small Small Medium Large
 X-Large 2X-Large 3X-Large 4X-Large

TM6 Hinge — Includes Extension Stops

Optional Flexion Stop Kit*

Correction *(Hinge Bar Will Be Pre-Bent)*

3° Standard Corrective Force
 5° Moderate Corrective Force
 7° Maximum Corrective Force

Replacement Sleeves *(Indicate Replacement Sleeve Size)*

X-Small Small Medium Large X-Large
 XX-Large XXX-Large XXXX-Large

Sizing Reference

	Knee Circumference	Calf Circumference	Thigh Circumference
X-Small	12" to 13"	12" to 13"	16" to 24"
Small	13" to 14"	13" to 15"	17" to 26"
Medium	14" to 15"	14" to 18"	18" to 28"
Large	15" to 17"	15" to 20"	20" to 29"
X-Large	17" to 19"	16" to 22"	22" to 31"
2X-Large	19" to 21"	18" to 24"	23" to 32"
3X-Large	21" to 22"	20" to 26"	25" to 34"
4X-Large	22" to 23"	21" to 27"	26" to 36"

For Multiple Brace (Stock Inventory) Orders

Quantity	Brace Size	
	X-Small	For stock braces, unless otherwise specified, the bars are pre-bent for 3 degrees of correction. Sleeves are universal (left or right leg). As such, a sleeve with a bar inserted in the pocket on the right side of the sleeve can be used as a right leg medial OA brace, or as a left leg lateral OA brace (with a medial hinge). A sleeve with a bar inserted in the pocket on the left side of the sleeve can be used as a left leg medial OA brace, or as a right leg lateral OA brace (with a medial hinge). If you need to use a brace for lateral OA, with a lateral hinge position, the bar must be re-bent. Instructions are provided with the brace to ensure the fitter has the hinge/bar assembly inserted into the appropriate side of the sleeve and bent in a direction that will provide the desired unloading forces. When ordering, you can request that we fill orders for a specific leg and compartment. But to minimize inventory, you can maintain one brace in each size and reposition (and potentially re-bend) the bar as needed.
	Small	
	Medium	
	Large	
	X-Large	
	2X-Large	
	3X-Large	
	4X-Large	
NOTE: Braces shipped for stock inventory will be packaged with the hinge inserted in the right pocket (right leg medial, left leg lateral).		

Please complete and fax this form to 800.798.2722 (24-hours a day). If you are calling in your order, this form indicates the options and information that will be required by our staff. For phone orders, please call 800.700.2722 between 6:00 a.m. and 4:00 p.m. (PST).

*Indicates additional charges apply

Thuasne USA

4615 Shepard St., Bakersfield, CA, 93313
 Phone: 800.432.3466 or 661.837.1795; Fax: 800.798.2722
 www.ThuasneUSA.com

